

CARE OF INSTRUMENTS

		Dilution Rate	Soak Time
DAILY USE:	Rapid	Light - 5mls	4 - 5 min
	Multi-Enzyme	Moderate - 7.5mls	8 min
		Heavy - 10mls	8 - 9 min
WEEKLY USE:	Surgislip	1 part to 10	
MONTHLY USE:	Surgistain	1 part to 7	

- Please read and follow all instructions on the cleaning products you are using.
- When new instruments are purchased wash and dry them also lubricate with either Aesculap Oil or Aerosol.



Shoof International Ltd, Private Bag 522, Cambridge 3450
Call free on 0800 800-801 • Fax free on 0800 264-625



310108