

Products manufactured in New Zealand by Shoof International Ltd. Index on back cover.



### Airline Sprayer Dairy

complete **214 675**

The Farmhand Airline Sprayer was invented and developed by a New Zealand dairy farmer. Now manufactured in New Zealand by Shoof International Ltd. This sprayer can be used in many industries for achieving perfectly atomised spray of all manner of liquid or cream products, even products of quite high viscosity. For the dairy industry the sprayer can be used for perfect application of many teat spray and teat barrier products that are unsprayable in any other teat spray system. Not only does the Farmhand Airline Sprayer give optimal atomisation and teat or udder coverage, but it generally also does this using less remedy. So it can easily recoup its cost when used with high-value teat barrier products. For the dairy industry the sprayer can be used for perfect application of many teat spray and teat barrier products that are unsprayable in any other teat spray system. Not only does the Farmhand Airline Sprayer give optimal atomisation and teat or udder coverage, but it generally also does this using less remedy. So it can easily recoup its cost when used with high-value teat barrier products. The wand tip incorporates a hook for quick and easy hang-ups between uses. The whole unit is very simple and convenient to operate in any dairy. The principal of operation is that the air flow through the handle draws product with it to be atomised as it passes through the nozzle. The spray pattern is a filled-in cone.



### Mastitis Test X-Spurt

black complete **220 910** **NEW**  
blue complete **204 475**  
white complete **210 553**

The X-Spurt paddle makes RMT testing a fast and efficient one-hand operation as the 300ml (10fl. oz) X-Spurt reagent reservoir is also the handle for the paddle. Milk is entered in the paddle chambers and the paddle tilted in the normal way to reduce the samples to the usual 3ml (0.1fl.oz). A quick squeeze of the X-Spurt handle then enters equal amounts of reagent at the same time to all four samples. Swirling the paddle mixes the samples for usual coagulation recognition. The bottle capacity is enough to RMT test about 15 cows. The X-Spurt is available in black, blue, or white.




Sheepmat on a plywood base in a sheep yard.

### Hoofmat™ Sheepmat

premium (blue) **217 592**

The Sheepmat is manufactured in the same tradition as the famous Hoofmat for cows. Premium model uses 31mm foam in a premium quality blue fabric outer. This should suit most existing race systems but the mat must be placed on concrete, wood, or other impervious surface. For successful treatment use in a situation where the animals will not be walking in dirty conditions, or for long distances, immediately after crossing the mat. Management systems may have to be developed to enable lame or at-risk animals to be treated to cross the mat at least once per day. Sheepmat should be used with about 10 litres of 2% copper sulphate solution, or other remedy as may be recommended by your veterinarian. Size of the sheepmat is 1.5m long x 0.46m wide.



### Quartermilkier 8 Litre

complete **207 556**  
repl't top assy cpt **207 564**  
repl't o-ring only **205 282**  
repl't rubber tube each **205 286**  
repl't connector **205 279**

For separating one quarter from the other three quarters in cases of colostrum, haemorrhage and high cell count (and antibiotic contamination where regulations allow).

Guaranteed unbreakable. Easy to clean. All antibiotic contaminated milk should be dumped, but the other three quarters should be suitable for calf rearing (or for sale where regulations allow). Our 8 litre Quartermilkier features a low-profile non-capsize design and specially moulded soft rubber tubes which allow the cow plenty of free movement whilst milking. Large 9cm neck opening. All parts available. Instructions for use included.



### Mastitis Test Super Paddle

black **211 537** **NEW**  
blue **204 473**

RMT paddle designed and manufactured by Shoof with special features over other paddles: smooth and comfortable handle, tilt lines for both left and right hand use, slightly domed cup-bottoms to enhance coagulation visibility, dark blue colour for maximum milk visibility. Quality long-life flexible polypropylene manufacture. Available in other colours if required, by indent order.