

BOOK....



"Goatkeeping 101"

210 615 \$85.95

As its name implies, this book has all the basics, plus huge amounts of advice and information of great value for the beginning breeder, and for the 'old hand'. Soft cover, 256 pages, with comprehensive index and lots of photos and diagrams.

PRICE REDUCED



"The Goatkeepers Veterinary Book"

200 854 ~~\$89.95~~ \$82.95

A standard veterinary reference work for all goat keepers. Peter Dunn (author) is a vet surgeon and goatkeeper/breeder in France. Includes a goatkeeper's calendar of particular interest in preventive medicine planning. Hardback, 210 pages, colour.

"Raising Milk Goats Successfully"

200 852 \$59.95



A good book for the beginner, covering basic management, nutrition, rations, housing, kid care, breeding, health, showing, plus sections on how to choose sires and dams. Lots of information. Paperback 172 pages.

"Goat Health Handbook"

200 845 \$89.95



If you raise goats you absolutely should have a copy of this book even if you have no others! Divided into 5 sections: diagnostic guides, disease description, therapy, treatment techniques, birth and newborns. Lots of tables and illustrations. Spiral bound paperback 123 pages.

"The Veterinary Book for Sheep Farmers"

200 857 \$115.00



The world's leading veterinary manual for sheep farmers. Produced in a style conducive to use by the average sheep breeder. Although produced in the UK and using that farming model, the information is directly applicable to sheep farming worldwide. Full of black and white photographs and detailed diagrams. Hard cover, 690 pages. Good index.



Goat Bell

65mm 203 227 \$21.95

A distinctive sound in the countryside in many northern hemisphere countries is the ubiquitous cow, goat and sheep bell. Now we are offering these here. Our first model is a simple brass coloured steel goat bell, 65mm wide, and 75mm deep. Suitable for a collar up to 40mm wide. Use one of our goat collars, or get a nice leather one to look the part. These bells have a distinctive sound.



PRICE REDUCED

Neck Chain Tags

blank x 10-pack 204 685 ~~\$26.95~~ \$22.95

numbered x10 204 686i ~~\$38.95~~ \$34.95

(prices on application for larger quantities)

Unbreakable solid nylon tags available blank (white), or numbered (red). The engraved tags are permanent numbers that can not wear or be chewed off. Hole is 19mm suitable for plastic chains, or kid and doe collars. Tag is a big 9cm x 5cm. Numbered tags are supplied on indent-order basis only.



Goat Collar Nylon

kid each (see chart) \$7.95

doe each (see chart) \$10.95

buck each (see chart) \$11.95

	Kid	Doe	Buck
black	203 236	black 203 232	black 203 228
blue	203 237	blue 203 233	blue 203 229
red	203 240	red 203 235	red 203 231
green	203 238	green 203 234	green 203 230

Long lasting nylon collars, flexible and reliable. Triple stitched to take hard wearing. Extra welded O-ring for tie-ups.

Kid: 16mm wide x 40cm long to suit up to 3-4 months.

Doe: 25mm wide x 55cm long with extra eyelets to suit smaller animals.

Buck: 25mm wide x 65cm long also fits large doe.



Sheepmat on a plywood base in a sheep yard.

Hoofmat™ Sheepmat

standard (black) 203 660 ~~\$145.00~~ \$99.95

premium (blue) 203 658 \$155.00

The Sheepmat is manufactured in the same tradition as the famous Hoofmat for cows. The Standard model uses 25mm foam with black outer polyfabric and the Premium model uses 31mm foam in a premium quality blue fabric outer. This should suit most existing race systems but the mat must be placed on concrete, wood, or other impervious surface. For successful treatment use in a situation where the animals will not be walking in dirty conditions, or for long distances, immediately after crossing the mat. Management systems may have to be developed to enable lame or at-risk animals to be treated to cross the mat at least once per day. Ideally the Sheepmat should be used in advance to harden the horn tissue of the hoof and thus negate the requirement for intensive treatment after deterioration. Even passing over the Sheepmat once or twice a week for several weeks before critical periods should be effective in reducing hoof disorders. Sheepmat should be used with about 10 litres of 20% copper sulphate solution, or other remedy as may be recommended by your veterinarian. Size of the sheepmat is 1.5m long x 0.46m wide.



Footrot Trough 2.4m

211 416 \$225.00

Designed to fit in a raceway for easy footrot treatment of sheep and goats. Sloping ends prevent leg damage. Ribbed base prevents slipping. Strong plastic manufacture will not rot or crack. Easily moved and emptied. Dimensions: 2.4m L x 34cm W x 20cm D. Capacity approx 70 litres.